

Chart a full test battery without sending patients to the cloud.

A self-hosted workstation for private-practice and school audiologists: patients, visits, visual audiograms, tympanometry with ipsi/contra reflexes, speech scores, and a clean printout. Runs on a Windows PC you already own.


Click-to-plot audiogram

Tymp + 500–4k reflexes

Speech / WRS

Print to PDF

No cloud PHI



Live product: air and bone, masked and unmasked, ASHA symbols, 3- and 4-frequency PTA. Click the chart or type the table — they stay in sync.

One visit. Audiogram, tymps, speech, print.

The screenshot shows a 'Today's visits' table with columns for patient name, date, and status. Below it is a 'Recent patients' table with columns for patient name, date, and status. The tables are populated with demo data.

Today's visits		
patient	date	status
John Doe	2024-01-15	Completed
Jane Smith	2024-01-15	In Progress
Bob Johnson	2024-01-15	Completed
Alice Brown	2024-01-15	In Progress
Charlie White	2024-01-15	Completed
Diana Green	2024-01-15	In Progress
Frank Black	2024-01-15	Completed
Grace King	2024-01-15	In Progress
Henry Lee	2024-01-15	Completed
Ivy Scott	2024-01-15	In Progress
Jack Hall	2024-01-15	Completed
Karen Adams	2024-01-15	In Progress
Leo Baker	2024-01-15	Completed
Mia Clark	2024-01-15	In Progress
Noah Evans	2024-01-15	Completed
Olivia Foster	2024-01-15	In Progress
Peter Garcia	2024-01-15	Completed
Quinn Hill	2024-01-15	In Progress
Samuel King	2024-01-15	Completed
Tina Lewis	2024-01-15	In Progress
Victor Miller	2024-01-15	Completed
Wendy Moore	2024-01-15	In Progress
Xavier Nelson	2024-01-15	Completed
Yara Ortiz	2024-01-15	In Progress
Zoe Parker	2024-01-15	Completed

Clinic board: today's visits, recent patients, one-click new chart.



Tympanometry with Jerger type, schematic curve, IPSI and CONTRA at 500, 1k, 2k, and 4k Hz.

In the test, not after it

- Patient search, MRN, referring provider
- Save-as-you-go on every battery
- Speech: SRT, SDT, WRS (NU-6 / CID W-22), MCL, UCL, binaural
- Print: name and date of birth on every page; empty sections are dropped so you don't print blank tymp or speech pages
- Clinician and admin logins on a shared clinic PC

Windows install in minutes

Unzip, open **installer\windows**, double-click **install.bat**, pick a mode. No Node knowledge required.

CLIENT

Booth PC. Opens the clinic in a browser. Points at your server. No local database.

SERVER

Shared host. Stores charts. Start / Stop shortcuts. Other booths connect as clients.

DUAL

One-room clinic or laptop: this PC is both the server and the workstation.

Data stays on the clinic PC or LAN you control. This is not a cloud EHR and not a HIPAA certification — keep the SQLite file on encrypted clinic storage and do not expose the server to the public internet. Screenshots use fake demo patients.